

Installation Illustrated

Aquaresp Quick Install Guide
by
Morten

Download

The screenshot shows the GitHub repository page for `bigb8/AquaResp`. The repository is public and has 2 branches and 3 tags. The `Code` dropdown menu is open, showing options to clone the repository using HTTPS, SSH, or GitHub CLI, or to download the ZIP file. The repository description is "Automating Aquatic Respirometry. Focused on automating intermittent respirometry experiments for aquatic animals." The repository is licensed under GPL-3.0. The latest release is "New oxygen sensor and controls" published 2 days ago. The repository is also compatible with the new firmware from PyroScience.

bigb8/AquaResp: Automating / X +

https://github.com/bigb8/AquaResp

Search or jump to... Pull requests Issues Marketplace Explore

bigb8 / AquaResp Public

Unwatch 2 Unstar 1 Fork 2

<> Code Issues Pull requests Discussions Actions Projects Wiki Security Insights Settings

master 2 branches 3 tags

Go to file Add file Code

bigb8 Update README.md

Aquarasp 3 Create 20211102.png

.gitattributes Create .gitattributes

.gitignore Update .gitignore

LICENSE Create LICENSE

README.md Update README.md

_config.yml Set theme jekyll-theme-midnight 2 years ago

Clone

HTTPS SSH GitHub CLI

https://github.com/bigb8/AquaResp.git

Use Git or checkout with SVN using the web URL.

Open with GitHub Desktop

Download ZIP

About

Automating Aquatic Respirometry. Focused on automating intermittent respirometry experiments for aquatic animals.

Readme

GPL-3.0 License

Releases 3

New oxygen sensor and controls Latest 2 days ago

+ 2 releases

Packages

No packages published

Publish your first package

Environments 1


github pages


AquaResp

Automating Aquatic Respirometry. Focused on automating intermittent respirometry experiments for aquatic animals.


Now compatible with the new firmware from PyroScience. See and use [OxyLog](#)

Right click zip file and choose
Extract all







Move or copy to desktop




To install Python, first download by 1.vbs



Allow running file






Python installer is downloaded when you see
“Installer.exe”



Run 2.bat

Run




Wait



RESTART YOUR PC


just to be sure that Windows
has registered Python

To check libraries- Run CheckLibraries.py - either by clicking or via CMD



Make links on desktop (optional, but recommended)





AquaResp v.3

Experiment name

Start Experiment

Animal settings

Equipment and experiment settings

Flush: 90s, Wait: 30s, Measure: 60s

Standard SMR (Flush -> Wait -> Measure)

Number of chambers: 4

Sensor: Presens Fibox 3

Verification of pO2 values:

Channel 1: 0

Channel 2: 0

Channel 3: 0

Channel 4: 0




Administrer Skrivebord

Genrejsværditager

Denne pc > Skrivebord

Søg i Skrivebord

	Navn	Ændringsdato	Type	Størrelse
Skrivebord	AquaResp 3 ASAP	04-11-2021 20:19	Genvej	2 KB
Overførsler	AquaResp Plots	04-11-2021 20:19	Genvej	2 KB
Dokumenter				



Open AquaResp.

Check that oxygen sensors are live.

And that the pump works

N.B.: Sensor Software needs to be configured and logging.

